

Patient Functions

- [How To: Enroll a Patient](#)
- [How to: Unenroll a Patient](#)
- [How To: Use Bypass](#)
- [How To: Suspend Tamper](#)
- [How To: Locate Patient](#)

How To: Enroll a Patient

Welcome, [admin@mychildcloud.com](#) | 10 | Show

Powered by McRoberts Technologies

Visit us to

Home Patients s Admin Help Logout

Click Patients Tab

TAG SUMMARY

Status	Count
InUse	6
Available	44
LowBattery	21
Bypass	1

Inventory Target: **N/A** Inventory Total: 72


SUMMARY

Status	Count
Active	6
Bypass	1

UNACCEPTED ALARMS

Tag near open door - 02/11/2022 20:11:32	Level 3
Tag near open door - 02/11/2022 20:09:32	Level 2
Tag near open door - 02/11/2022 18:44:03	Level 1
Tag near open door - 02/09/2022 21:04:40	
Tag near open door - 02/09/2022 21:04:35	
Tag near open door - 02/09/2022 21:04:31	
Tag near open door - 02/09/2022 21:04:26	
Tag near open door - 02/09/2022 21:04:21	
Tag near open door - 02/09/2022 21:04:18	

dev.mychildcloud.com:8080/Patients/Directory/



The enrollment window will pop up. Select the tag you wish to use from the drop-down list.

Enroll Patient [X]

First Name

Middle Name

Last Name

Room

Gender

☐ Flight Risk

Unit

Tag

Floor

Enroll

© 2024 - McRoberts Technologies | build version: 1.0.0

Home Patients Tags Floorplans Reports Admin Help My Profile

Enter the patient's information. If you leave any required fields blank you will be prompted to fill them out. Only the "Middle Name" field is optional and can be left blank.

Enroll Patient

First Name

Middle Name

Last Name

Room

Unit

Tag

Floor

Sex ☒ Male ☐ Female

☒ Flight Risk

☒ Mother Required

☐ Court Order

☐ Other

Enroll



Enter Patient Information


Then Click Enroll

© 2024 - McRoberts Technologies | build version: 1.0.0

Patient will now have a status of "Active" and the total number of patients will be displayed on the dashboard. Other patient status are "Bypass" and "Tamper Suspended."

How to: Unenroll a Patient






After clicking "OK" you will receive a notification that the patient has been unenrolled. The page will then refresh and the patient will be removed from the patient directory. The tag status will change to "available" meaning it can be used on a new patient.

How To: Use Bypass


MyChild provides a ByPass that allows you to take a patient out of a protected area without setting off alarms. You must select a departing unit and a destination unit. These can be the same. The patient status changes from “Active” to “Bypass” when a patient is in ByPass.



Patient Directory

Enroll Patient


Ans: Or 47 47 X Apply




Patient Directory

[Enroll Patient](#)

Avail On X



Enter a reason for the Bypass and a time limit.



The screenshot shows a 'Patient Bypass' form with the following fields and annotations:

- Departing Unit:** Labor and Delivery (dropdown menu). An arrow points to this field from a red box containing the text: "Destination unit may be the same as the departing unit OR a different protected unit."
- Destination Unit:** Mother / Baby (dropdown menu).
- Time Required (in minutes):** 5 (input field with up/down arrows). An arrow points to this field from a red box containing the text: "Enter in the time you will be out of protected units."
- Reason for Bypass:** Moving (dropdown menu).
- Transporter Name:** Jill Smith (Optional) (text input field).
- Transporter Phone:** (248) 555-9999 (Optional) (text input field).
- Put Patient In Bypass** (button).

Once the information is complete, click the "Put Patient In Bypass" button. This will change the patient status to Bypass. You will also see this change on the dashboard under the Patient Summary.

A level one notification will appear on the upper right corner of the screen telling you the Bypass has started once the patient leaves through the first bypass designed door.

MyChild

MyChild

Home

Patients

Tags

Floorplans ▾

Reports ▾

Admin ▾

Help ▾

Logout

LABOR AND DELIVERY

McRoberts HOSPITAL

BYPASS STARTED

Patient Name:

Davis, Brad

Tag Number:

B2DE55

Location Start:

Labor and Delivery

Location End:

Labor and Delivery

Transporter name:

Nurse Jackie

Transporter phone #:

800-555-1234

Time of Alarm:

05/03/2022 17:43:09

Accept

Edit Bypass

Labor and Delivery

© 2022 - McRoberts Technologies | build version: 8122

When the time allotted for the bypass is almost expired, a “BYPASS EXPIRING” alarm is triggered.

Home Patients Tags Floorplans ▼ Reports ▼ Admin ▼ Help ▼ Logout

LABOR AND DELIVERY

McRoberts HOSPITAL

BYPASS EXPIRING

170 secs left


Patient Name: Davis, Brad
Tag Number: B2DE55
Location Start: Labor and Delivery
Location End: Labor and Delivery
Transporter name: Nurse Jackie
Transporter phone #: 800-555-1234
Time of Alarm: 05/03/2022 17:58:06

Accept **Edit Bypass**

Click "Accept" to clear the alarm from the screen. The timer will continue to run out until the bypass expires or is complete.

Click "Edit Bypass" to add more time to the bypass or edit the bypass details.

When the timer reaches zero, a "BYPASS EXPIRED" alarm is triggered. All bypass alarms (except BYPASS COMPLETE) allow you to edit the bypass. Changes can also be made on the patient card.




First, click the "Accept" button, then enter a message and click the green arrow to accept the alarm. The Bypass Expired alarm only happens if the time limit has run out and the patient has not passed through the second bypass designated door. If the patient arrives on time, you will receive a Bypass Complete alarm. This will show up in the upper right corner as a notification as it is a level one alarm. Click on the notification to view the alarm detail.



LABOR AND DELIVERY

McRoberts
HOSPITAL



Labor and Delivery

BYPASS COMPLETE

Patient Name:	Davis, Brad
Tag Number:	B2DE55
Location Start:	Labor and Delivery
Location End:	Labor and Delivery
Transporter name:	Nurse Jackie
Transporter phone #:	800-555-1234
Time of Alarm:	05/03/2022 17:44:29

Accept

When you complete a bypass by clicking “Complete,” you are acknowledging that the patient has arrived at the destination unit. The patient status will return to "Active".

For additional help, call 1-800-776-7328 and press 2 for technical support.

[Click here to download PDF.](#)

How To: Suspend Tamper

MyChild provides a suspend tamper function which allows you to remove a tag temporarily without setting off alarms. To remove a cut-band or sensing tag temporarily, put the patient in SUSPEND TAMPER. The patient status changes from “Active” to “Tamper Suspended.”

The screenshot displays the MyChild Patient Directory interface. At the top is a teal header with the MyChild logo and a navigation menu: Home, Patients, Tags, Floorplans, Reports, Admin, Help, and Logout. Below the header is the 'Patient Directory' section. It features a grid of patient cards. The first card on the left is a teal card with a white circle containing a plus sign and the word 'NEW'. Below this is a larger teal card with the MyChild logo and the text 'SECURITY FIRST'. To the right are two more teal cards. The first of these is for 'Davis, Brad', showing 'SEX: MALE', 'UNIT: PEDS', and 'TAG: D1F5A2'. The second is for 'Olt, Anna', showing 'SEX: FEMALE', 'UNIT: PEDS', and 'TAG: E150CD'. Both of these cards have a pink 'FUNCTIONS' button with a refresh icon. Below the 'Olt, Anna' card is a yellow bar with the text 'FLIGHT RISK'. At the bottom, there is a partial view of a card for 'Davis, Sarah' and another card partially visible on the right.

Patient Name	Sex	Unit	Tag	Status
Davis, Brad	MALE	PEDS	D1F5A2	Active
Olt, Anna	FEMALE	PEDS	E150CD	FLIGHT RISK
Davis, Sarah	FEMALE			

Complete Suspend Tamper details:

The image shows a mobile application interface. On the left, a 'Patient Suspend Tamper' dialog box is open over a patient card. The dialog box contains the following fields: 'Time Required (in minutes):' with a numeric input field, 'Reason for Suspend:' with a text input field, and 'Nurse's Name:' with a text input field. A 'Create' button is at the bottom of the dialog. The patient card in the background shows 'FEMALE PEDS' and 'F1D71B'. At the bottom of the card are two buttons: 'FUNCTIONS' with a refresh icon and 'FLIGHT RISK'. On the right, another patient card for 'Woods, Gloria' is shown. It has a 'NO PHOTO' placeholder and four buttons: 'BYPASS', 'SUSPEND TAMPER' (highlighted with a yellow box), 'UNENROLL', and 'LOCATE'.

Once the information is complete, click the Create button.
This will change the patient status to “Tamper Suspended.”

When the time allotted for the suspend tamper is almost expired, a “SUSPEND TAMPER EXPIRING” alarm is triggered. When the timer reaches zero, a “SUSPEND TAMPER EXPIRED” alarm is triggered. All suspend tamper alarms allow you to edit the suspend tamper. Changes can also be made on the patient card.

When you click “Complete” to complete a Suspend Tamper, you are acknowledging that the tag is back on the patient.

How To: Locate Patient