


Patient Functions

- [How To: Enroll a Patient](#)
- [How to: Unenroll a Patient](#)
- [What is Bypass?](#)
- [How To: Suspend Tamper](#)
- [How To: Locate Patient](#)

How To: Enroll a Patient



The screenshot shows the MyChildCloud web application interface. At the top is a teal header with the MyChildCloud logo. Below the header is a navigation bar with tabs: Home, Patients, Tags, Admin, Help, and Logout. A red callout box with the text "Click Patients Tab" and a red arrow points to the Patients tab. The main content area is divided into three sections: TAG SUMMARY, TAG SUMMARY, and UNACCEPTED ALARMS. The first TAG SUMMARY section contains a table with Status and Count columns, showing InUse (6), Available (44), LowBattery (21), and Bypass (1). Below this table are fields for Inventory Target (N/A) and Inventory Total (72). The second TAG SUMMARY section contains a table with Status and Count columns, showing Active (6) and Bypass (1). The UNACCEPTED ALARMS section contains a list of alarms, each with the text "Tag near open de".

MyChildCloud

Home Patients Tags Admin Help Logout

Click Patients Tab

TAG SUMMARY

Status	Count
InUse	6
Available	44
LowBattery	21
Bypass	1

Inventory Target: **N/A** Inventory Total: 72

TAG SUMMARY

Status	Count
Active	6
Bypass	1

UNACCEPTED ALARMS

- Tag near open de
- Tag near open de
- Tag near open de
- Tag near open de
- Tag near open de
- Tag near open de
- Tag near open de
- Tag near open de
- Tag near open de
- Tag near open de

dev.mychildcloud.com:8080/Patients/Directory/

Patient Directory



To add a
new patient,
click here



The screenshot shows the 'MyChild' Patient Directory interface. At the top, there is a navigation bar with links: Home, Patients, Tags, Floorplans, Reports, Admin, Help, and Logout. The main section is titled 'Patient Directory' and displays a grid of patient cards. One card is highlighted with a '+ NEW' button. An 'Enroll Patient' modal is open, showing fields for First Name, Middle Name, Last Name, Department, Tag, Floor, Sex (Male/Female), and a checkbox for Flight Risk. A red box highlights the 'Enroll' button at the bottom of the modal. A red arrow points to the 'Enroll' button from the right side of the image.

MyChild

Home Patients Tags Floorplans Reports Admin Help Logout

Patient Directory

Enroll Patient

First Name

Middle Name

Last Name

Department Please select a unit

Tag Please select a tag

Floor Please select a floor

Sex ☒ Male ☐ Female

☐ Flight Risk

Enroll

Patient will now have a status of "Active" and the total number of patients will be displayed on the dashboard. Other patient status are "Bypass" and "Tamper Suspended."

How to: Unenroll a Patient

The screenshot shows the MyChild Patient Directory interface. The header includes the MyChild logo and a navigation bar with links: Home, Patients, Tags, Floorplans, Reports, Admin, Help, and Logout. The main content area is titled "Patient Directory" and displays a grid of patient cards. The first card is a "NEW" button. The other cards show patient details: Davis, Brad (MALE, DefaultUnit, D1F5A2), Olt, Anna (FEMALE), Davis, Sarah (FEMALE), and Tag, Test (MALE). A red callout box highlights the "FUNCTIONS" button on the Davis, Brad card, with the text: "Access all patient functions by clicking the pink functions button on the card".

MyChild

Home Patients Tags Floorplans Reports Admin Help Logout

Patient Directory

NEW

Davis, Brad

SEX: MALE
UNIT: DefaultUnit
TAG: D1F5A2

FUNCTIONS

Olt, Anna

SEX: FEMALE
UNIT:
TAG:

Davis, Sarah


SEX: FEMALE
UNIT:
TAG:

Tag, Test

SEX: MALE
UNIT:
TAG:

Access all patient functions by clicking the pink functions button on the card

Patient Directory



**Click
UNENROLL
button.**



Image not found or type unknown



After clicking "OK" you will receive a notification that the patient has been unenrolled. The page will then refresh and the patient will be removed from the patient directory.

What is Bypass?

MyChild provides a ByPass Mode that allows you to take a patient out of a protected area without setting off alarms. You must select a departing unit and a destination unit. These can be the same. The patient status changes from “Active” to “Bypass” when a patient is in ByPass Mode.

The Bypass Process

Put patient into Bypass Mode and select how much time you want the Bypass Mode to be in effect.



Bypass time begins when you pass the **first bypass-designated door**



Bypass Mode ends when you pass the **first bypass-designated door a second time.**





MyChild

[Home](#)[Patients](#)[s ▼](#)[Admin ▼](#)[Help ▼](#)[Logout](#)

Click
Patients
Tab

TAG SUMMARY

Status	Count
InUse	6
Available	44
LowBattery	21
Bypass	1

Inventory Target: **N/A** Inventory Total: 72

SUMMARY

Status	Count
Active	6
Bypass	1

UNACCEPTED ALARMS

Tag near open d

Tag near open d

Tag near open d

Tag near open d

Tag near open d

Tag near open d

Tag near open d

Tag near open d

Tag near open d




Patient Directory



Access all patient functions by clicking the pink functions button on the patient card.


Patient Directory



**Click Bypass
to continue**



Enter a reason for the Bypass, as well as your name



The image displays two components of a patient bypass system. On the left is a patient card for 'Davis, Brad' with a 'NO PHOTO' placeholder and a refresh icon. It features four teal buttons: 'BYPASS' (highlighted with a yellow border), 'SUSPEND TAMPER', 'UNENROLL', and 'LOCATE'. On the right is the 'Patient Bypass' form, which includes fields for 'Departing Unit' (set to 'PEDS'), 'Destination Unit' (a dropdown menu), 'Time Required (in minutes)' (a numeric input field with a red arrow pointing to its spinner controls), 'Reason for Bypass', 'Transporter Name', and 'Transporter Phone'. A 'Create' button is located at the bottom right of the form.

Once the information is complete, click the Create button.
This will change the patient status to “Bypass.”

A level one notification will appear on the upper right corner of the screen telling you the Bypass has started once the patient leaves though the first bypass designed door.

MyChild

MyChild

Home

Patients

Tags

Floorplans ▾

Reports ▾


Admin ▾

Help ▾

Logout

LABOR AND DELIVERY

McRoberts HOSPITAL



Labor and Delivery

BYPASS STARTED

Patient Name:

Davis, Brad

Tag Number:

B2DE55

Location Start:

Labor and Delivery

Location End:

Labor and Delivery

Transporter name:

Nurse Jackie

Transporter phone #:

800-555-1234

Time of Alarm:

05/03/2022 17:43:09

Accept



Edit Bypass

© 2022 - McRoberts Technologies | build version: 8122

When the time allotted for the bypass is almost expired, a “BYPASS EXPIRING” alarm is triggered.

Home Patients Tags Floorplans ▼ Reports ▼ Admin ▼ Help ▼ Logout

LABOR AND DELIVERY



Labor and Delivery

BYPASS EXPIRING

170 secs left


Patient Name: Davis, Brad
Tag Number: B2DE55
Location Start: Labor and Delivery
Location End: Labor and Delivery
Transporter name: Nurse Jackie
Transporter phone #: 800-555-1234
Time of Alarm: 05/03/2022 17:58:06

[Accept](#) [Edit Bypass](#)

Click “Accept” to clear the alarm from the screen. The timer will continue to run out until the bypass expires or is complete.

Click “Edit Bypass” to add more time to the bypass or edit the bypass details.

When the timer reaches zero, a “BYPASS EXPIRED” alarm is triggered. All bypass alarms (except BYPASS COMPLETE) allow you to edit the bypass. Changes can also be made on the patient card.




First, click the "Accept" button, then enter a message and click the green arrow to accept the alarm. The Bypass Expired alarm only happens if the time limit has run out and the patient has not passed through the second bypass designated door. If the patient arrives on time, you will receive a Bypass Complete alarm. This will show up in the upper right corner as a notification as it is a level one alarm. Click on the notification to view the alarm detail.



LABOR AND DELIVERY

McRoberts
HOSPITAL



BYPASS COMPLETE

Patient Name: Davis, Brad
Tag Number: B2DE55
Location Start: Labor and Delivery
Location End: Labor and Delivery
Transporter name: Nurse Jackie
Transporter phone #: 800-555-1234
Time of Alarm: 05/03/2022 17:44:29

Accept

Labor and Delivery

When you complete a bypass by clicking “Complete,” you are acknowledging that the patient has arrived at the destination unit. The patient status will return to "Active".

For additional help, call 1-800-776-7328 and press 2 for technical support.

[Click here to download PDF.](#)

How To: Suspend Tamper

MyChild provides a suspend tamper function which allows you to remove a tag temporarily without setting off alarms. To remove a cut-band or sensing tag temporarily, put the patient in SUSPEND TAMPER. The patient status changes from “Active” to “Tamper Suspended.”

The screenshot displays the MyChild Patient Directory interface. At the top is a teal header with the MyChild logo and a navigation menu: Home, Patients, Tags, Floorplans, Reports, Admin, Help, and Logout. Below the header is the 'Patient Directory' section. It features a grid of patient cards. The first card on the left is a teal card with a white circle containing a plus sign and the word 'NEW'. Below this is a white box with the MyChild logo and the text 'SECURITY FIRST'. To the right are two more teal cards. The first of these is for 'Davis, Brad', showing 'SEX: MALE', 'UNIT: PEDS', and 'TAG: D1F5A2'. The second is for 'Olt, Anna', showing 'SEX: FEMALE', 'UNIT: PEDS', and 'TAG: E150CD'. Both have a pink 'FUNCTIONS' button with a refresh icon. Below the Anna Olt card is a yellow bar with the text 'FLIGHT RISK'. At the bottom, the top of a card for 'Davis, Sarah' is visible, showing 'SEX: FEMALE'.

Complete Suspend Tamper details:

The image shows a mobile application interface. On the left, a 'Patient Suspend Tamper' dialog box is open over a patient card. The dialog box contains the following fields: 'Time Required (in minutes):' with a numeric input field, 'Reason for Suspend:' with a text input field, and 'Nurse's Name:' with a text input field. A 'Create' button is at the bottom of the dialog. The patient card in the background shows 'FEMALE PEDS' and 'F1D71B'. At the bottom of the card are two buttons: 'FUNCTIONS' with a refresh icon and 'FLIGHT RISK' in a yellow bar. On the right, another patient card for 'Woods, Gloria' is shown. It has a 'NO PHOTO' placeholder and four buttons: 'BYPASS', 'SUSPEND TAMPER' (highlighted with a yellow box), 'UNENROLL', and 'LOCATE'.

Once the information is complete, click the Create button.
This will change the patient status to “Tamper Suspended.”

When the time allotted for the suspend tamper is almost expired, a “SUSPEND TAMPER EXPIRING” alarm is triggered. When the timer reaches zero, a “SUSPEND TAMPER EXPIRED” alarm is triggered. All suspend tamper alarms allow you to edit the suspend tamper. Changes can also be made on the patient card.

When you click “Complete” to complete a Suspend Tamper, you are acknowledging that the tag is back on the patient.

How To: Locate Patient