

MyChild6-Clinical

- [How To Apply a Tag](#)
 - [Cut-Band Tag Application, Removal, and Disinfection Guide](#)
 - [Bracelet Tag Application, Removal, and Disinfection Guide](#)
 - [Tummy Tag Application, Removal, and Disinfection Guide](#)
- [Patient Functions](#)
 - [How To: Enroll a Patient](#)
 - [How to: Unenroll a Patient](#)
 - [How To: Use Bypass](#)
 - [How To: Suspend Tamper](#)
 - [How To: Locate Patient](#)
- [Tag Management](#)
 - [How To: Disinfect Tags](#)
 - [STEP Program](#)
 - [How To: Store and Test Tags](#)
 - [Unable to Assign a Tag to a Patient - Troubleshooting](#)
 - [How to: Store Radio Frequency Tags](#)
 - [Tag Warranty/Buy-Back](#)
- [General Functions](#)
 - [How To: Change your Password](#)
 - [How To: Enable 2FA](#)
 - [How To: Update SMS settings](#)
- [Alarms](#)
 - [Alarm Overview](#)

- Patient Alarms
 - Device Alarms
 - System Alarms
 - Alarm Chart
- Spec Sheets
 - Cut-Band Spec Sheet
 - Tummy Tag Spec Sheet
 - Bracelet Tag Spec Sheet
 - Wireless Receiver Spec Sheet
 - Bypass Keypad Spec Sheet
- STEP Customers - Instructions
- Changelogs
 - v6.0.7-SP2
 - 6.0.7-SP1

How To Apply a Tag

Cut-Band Tag Application, Removal, and Disinfection Guide



Cut-E Remo

App

1. **ASSIGN** the tag to a patient i
(See How to Store and Test a MyChild
2. **THREAD** the cut-band UP thr
of tag (white side with writin



3. **SLIP** the loop over baby's foot



Cut-Band Removal

Removing

1. **DISCHARGE** the tag in the MyChild scanner.
2. **CUT** the band as shown. Position the tag in the scanner.
3. **CLEAN** the tag thoroughly in accordance with the manufacturer's instructions.
4. **PLACE** tag in foil pouch for storage.

2



[Click here to download the Cut Band Application Guide!](#)

Bracelet Tag Application, Removal, and Disinfection Guide

Bracelet Tag Application, Removal, and Disinfection Guide

Note: Applying a Bracelet Tag to a patient is done **after** the Tag's battery is tested and **after** the Tag has been assigned to the patient.

For more information on Tag battery testing, see *How to Store and Test a MyChild Security Tag*. For more information on assigning a Tag to a patient, see *Bracelet Tag Comprehensive Guide PowerPoint*.

Bracelet Tag Application

1. Feed the band through the two slots on the top of the tag (side with writing).

You may use any number of common bands.

2. Fasten the band around the limb to form a loop.
3. The Tag should be "just snug." Verify that the band and Tag do not restrict the patient's circulation or cause chafing.
4. Trim the excess band material if desired.
5. The Bracelet Tag can face either way. The benefit to having the top face up is you can view the Tag number.
6. Use Surgilast, if desired, to keep Tag in place.

Bracelet Tag Removal

Note: Remove the Tag **after** discharging the patient in the MyChild software.

To remove the tag, cut the band with scissor and slide the band out of the tag.

Bracelet Tag Disinfection Recommendations

- Thoroughly clean Tag with disinfectant wipe. McRoberts uses SaniZide Pro 1 disinfectant wipes which are broad-spectrum, EPA-registered, hospital-grade surface disinfecting wipes that kill the Novel Coronavirus (COVID-19) and meet OSHA (Occupational Safety and Health Administration) requirements regarding bloodborne pathogens and personal protection.
- Use a small soft-bristled brush to scrub each tag with germicidal solution until there is no visible evidence of blood, body fluids, or material on the Tag. McRoberts uses SaniZide Pro 1, which is a broad-spectrum, EPA-registered, hospital-grade surface disinfectant that kills the Novel Coronavirus (COVID-19) and meets OSHA (Occupational Safety and Health Administration) requirements regarding bloodborne pathogens and personal protection.
- Spray each Tag continuously with germicidal solution to ensure a minimum of 1 minute “wet time” with the solution.
- Wipe dry with a disposable, lint-free towelette.
- Spray Tags with alcohol (70%) to remove any chemical residue.
- Wipe dry with a disposable, lint-free towelette. McRoberts uses Wypall towelettes.

How To Apply a Tag

Tummy Tag Application, Removal, and Disinfection Guide



Applying a Tummy Tag to a
tested and *after* the Tag has
information on Tag battery
Security Tag. For more infor
Tummy Tag Comprehensive G

1. POSITION THE UMBILICA

- The MyChild Tummy
on umbilical cord cle



Tummy Tag and Disinfection

Tag Disinfection

- Thoroughly clean tag with disinfectant which are broad-spectrum, EPA-registered, effective against Novel Coronavirus (COVID-19) and meet requirements regarding bloodborne pathogens.
- Use a small soft-bristled brush to scrub evidence of blood, body fluids, or material with broad-spectrum, EPA-registered, hospital-grade disinfectant (COVID-19) and meets OSHA (Occupational Safety and Health) requirements regarding bloodborne pathogens and personal protective equipment.
- Spray each tag continuously with germicide solution with the solution.
- Wipe dry with a disposable, lint-free towel.
- Spray tags with alcohol (70%) to remove any remaining disinfectant.
- Wipe dry with a disposable, lint-free towel.

[Click here to download the Tummy Tag Application, Removal, and Disinfection Guide!](#)

Patient Functions

How To: Enroll a Patient

The screenshot shows the MyChildCloud web application interface. At the top, there is a header bar with the MyChildCloud logo on the left and user information on the right, including a welcome message for 'admin@mychildcloud.com', a notification bell with '10' alerts, and a 'Logout' button. Below the header is a navigation menu with tabs: 'Home', 'Patients', 'Admin', 'Help', and 'Logout'. The 'Patients' tab is highlighted with a red box and a red arrow pointing to it, with a text box saying 'Click Patients Tab'. The main content area is divided into three sections: 'TAG SUMMARY', 'SUMMARY', and 'UNACCEPTED ALARMS'. The 'TAG SUMMARY' section contains a table with columns 'Status' and 'Count', showing 'InUse' (6), 'Available' (44), 'LowBattery' (21), and 'Bypass' (1). Below this table, it shows 'Inventory Target: N/A' and 'Inventory Total: 72'. The 'SUMMARY' section contains a similar table with 'Status' and 'Count', showing 'Active' (6) and 'Bypass' (1). The 'UNACCEPTED ALARMS' section contains a list of alarms, each with a timestamp and a location, such as 'Tag near open door - 02/11/2022 20:11:32'. To the right of the alarms, there are three colored buttons: 'Level 3' (red), 'Level 2' (orange), and 'Level 1' (yellow). The URL 'dev.mychildcloud.com:8080/Patients/Directory/' is visible at the bottom left.

Welcome, [admin@mychildcloud.com](#) | 🔔 10 | [Show](#)

Powered by McRoberts Technologies  [Visit us to](#)

Home Patients **Click Patients Tab** Admin Help Logout

TAG SUMMARY

Status	Count
InUse	6
Available	44
LowBattery	21
Bypass	1

Inventory Target: **N/A** Inventory Total: 72


SUMMARY

Status	Count
Active	6
Bypass	1

UNACCEPTED ALARMS

Tag near open door - 02/11/2022 20:11:32	Level 3
Tag near open door - 02/11/2022 20:09:32	Level 2
Tag near open door - 02/11/2022 18:44:03	Level 1
Tag near open door - 02/09/2022 21:04:40	
Tag near open door - 02/09/2022 21:04:35	
Tag near open door - 02/09/2022 21:04:31	
Tag near open door - 02/09/2022 21:04:26	
Tag near open door - 02/09/2022 21:04:21	
Tag near open door - 02/09/2022 21:04:18	

dev.mychildcloud.com:8080/Patients/Directory/



The enrollment window will pop up. Select the tag you wish to use from the drop-down list.

Enroll Patient [X]

First Name

Middle Name

Last Name

Room

Gender

☐ Flight Risk

Unit

Tag

Floor

© 2024 - McRoberts Technologies | build version: 1.0.0

Home Patients Tags Floorplans Reports Admin Help My Profile

Enter the patient's information. If you leave any required fields blank you will be prompted to fill them out. Only the "Middle Name" field is optional and can be left blank.

Enter Patient Information

Then Click Enroll

Enroll Patient

First Name

Middle Name

Last Name

Room

Unit

Tag

Floor

Sex ☒ Male ☐ Female

☒ Flight Risk

☒ Mother Required



☐ Court Order


☐ Other

© 2024 - McRoberts Technologies | build version: 1.0.0

Patient will now have a status of "Active" and the total number of patients will be displayed on the dashboard. Other patient status are "Bypass" and "Tamper Suspended."

How to: Unenroll a Patient






After clicking "OK" you will receive a notification that the patient has been unenrolled. The page will then refresh and the patient will be removed from the patient directory. The tag status will change to "available" meaning it can be used on a new patient.

How To: Use Bypass


MyChild provides a ByPass that allows you to take a patient out of a protected area without setting off alarms. You must select a departing unit and a destination unit. These can be the same. The patient status changes from “Active” to “Bypass” when a patient is in ByPass.



Patient Directory

Enroll Patient

Ans: Or 47 47 X Apply




Access all patient functions by clicking the pink **FUNCTIONS** button on the card.

Patient Directory


[Enroll Patient](#)

Avail On X



Click Bypass
to continue

Enter a reason for the Bypass and a time limit.



The form is titled "Patient Bypass" and contains the following fields:

- Departing Unit: Labor and Delivery
- Destination Unit: Mother / Baby
- Time Required (in minutes): 5
- Reason for Bypass: Moving
- Transporter Name: Jill Smith (Optional)
- Transporter Phone: (248) 555-9999 (Optional)

At the bottom is a button labeled "Put Patient In Bypass".


Destination unit may be the same as the departing unit OR a different protected unit.

Enter in the time you will be out of protected units.

Bypass Mode starts when the tag passes a bypass-designated door in the departing unit. Bypass Mode ends when the tag in bypass passes through a designated door in the destination unit.

Once the information is complete, click the "Put Patient In Bypass" button. This will change the patient status to Bypass. You will also see this change on the dashboard under the Patient Summary.

A level one notification will appear on the upper right corner of the screen telling you the Bypass has started once the patient leaves through the first bypass designed door.



The dashboard shows a map of the Labor and Delivery - DefaultFloor. A yellow box labeled "BYPASS STARTED" is overlaid on the map, containing the following information:

- Full Name: Smith, Mike
- Tag Number: E1F418
- Assigned Room: 201
- Location: West Door
- Time Of Alarm: 02/27/2025 09:20:10 AM

Below the map, there is a "Bypass Started" notification box with the following details:

- Name: Rhodes, BG Casey
- Assigned Room: Not Set
- Tag Number: 82F866
- Alarm Time: 5/15/2024 2:37:49 PM

A red arrow points from the "View Alarm" button in the notification box to the "BYPASS STARTED" box on the map.

This alarm indicates the patient has gone through the first door and the countdown has begun.

When the time allotted for the bypass is almost expired, a “BYPASS EXPIRING” alarm is triggered.

McRoberts HOSPITAL

LABOR AND DELIVERY

BYPASS EXPIRING

170 secs left


Patient Name: Davis, Brad
Tag Number: B2DE55
Location Start: Labor and Delivery
Location End: Labor and Delivery
Transporter name: Nurse Jackie
Transporter phone #: 800-555-1234
Time of Alarm: 05/03/2022 17:58:06

Accept **Edit Bypass**

Click “Accept” to clear the alarm from the screen. The timer will continue to run out until the bypass expires or is complete.

Click “Edit Bypass” to add more time to the bypass or edit the bypass details.

When the timer reaches zero, a “BYPASS EXPIRED” alarm is triggered. All bypass alarms (except BYPASS COMPLETE) allow you to edit the bypass. Changes can also be made on the patient card.




First, click the "Accept" button, then enter a message and click the green arrow to accept the alarm. The Bypass Expired alarm only happens if the time limit has run out and the patient has not passed through the second bypass designated door. If the patient arrives on time, you will receive a Bypass Complete alarm. This will show up in the upper right corner as a notification as it is a level one alarm. Click on the notification to view the alarm detail.



LABOR AND DELIVERY

McRoberts
HOSPITAL



BYPASS COMPLETE

Patient Name: Davis, Brad
Tag Number: B2DE55
Location Start: Labor and Delivery
Location End: Labor and Delivery
Transporter name: Nurse Jackie
Transporter phone #: 800-555-1234
Time of Alarm: 05/03/2022 17:44:29

Accept

Labor and Delivery


When you complete a bypass by clicking “Complete,” you are acknowledging that the patient has arrived at the destination unit. The patient status will return to "Active".

For additional help, call 1-800-776-7328 and press 2 for technical support.

[Click here to download PDF.](#)

How To: Suspend Tamper

Suspend Tamper Alarms: MyChild p
remove a cut-band or sensing tag tem
into SUSPEND TAMPER mode, the pat



When the time allotted for the suspend EXPIRING” alarm is triggered. When it is triggered. All suspend tamper alarm be made on the patient card.

When you click Complete to complete back on the patient.

Alarm Name	Level	Alarm Description
SUSPEND TAMPER STARTED	1	This alarm notifies a patient has been Suspend Tamper that the Suspend timer has started The patient status from “Active” to “Suspended.”
SUSPEND	1	This alarm appears

When you click “Complete” to complete a Suspend Tamper, you are acknowledging that the tag is back on the patient.

How To: Locate Patient

Tag Management

How To: Disinfect Tags



Tag Disinfecting Procedure

Step Number	Procedure
Step One	Separate the tags into batches of no more than 10 tags per tray to increase surface area for spray solution to ensure each tag is fully wet. The tags must have been fully coated with germicidal solution. <i>McRoberts uses SaniZide Pro 1, which is a hospital-grade surface disinfectant that is EPA registered and meets OSHA (Occupational Safety and Health) requirements regarding bloodborne pathogens.</i>
Step Two	After the tags have been wet for one minute, use a brush to scrub each tag until there is no visible debris material on the tag. <i>If cleaning a Cut-Band Tag, use caution. Do not scrub the tag. Gently clean the teeth on a Cut-Band Tag. Clean the bottom of teeth. All debris, such as blood, must be removed.</i>

STEP Program

STEP - Simple Tag Exchange Program

Never buy a tag! Let McRoberts manage your tags for you.

STEP is a program in which McRoberts manages your tag inventory, testing and disinfecting and provides tags for per-patient use with the MyChild infant protection system. You use tags once then return them to McRoberts for testing and disinfection on a weekly basis. STEP integrates seamlessly with clinical practices and relieves hospitals from the burden of cleaning tags, managing inventory, and the purchasing process when tags expire.

The birth census of each hospital determines the number of tags that are assigned to a hospital's Total Tag Population. The initial shipment of tags to the hospital contains a three-week supply. Each week the hospital sends its dirty tags to McRoberts and McRoberts ships back the same number of tested and disinfected tags. The program is designed so that at any time, a third of the tags are in use, a third are in the hospital's inventory and a third are in transit or in cleaning.

The foil pouches that the tags come in should be kept and used to send dirty tags back to McRoberts in order to preserve their battery life and get them off the protected unit without alarming.

Already own your tags? McRoberts will give a \$30 per-tag credit for tags that are not dead or have a low battery. McRoberts will also always take back any dead tag for proper disposal.

Tags are received, tested and disinfected at McRoberts' Fulfillment Center at:

McRoberts Technologies
409 Joyce Kilmer Avenue
Suite 111
New Brunswick, NJ 08901

The STEP contract is for twelve months and invoicing is monthly, based on anticipated birth census. Clients who opt to pay for the entire twelve months in advance enjoy a ten percent discount.

[Click Here](#) to download the PDF

For questions about the STEP Program, contact Tag Fulfillment at 800.776.7328, Extension 2104 or Fulfillment@McRobertsTech.Com.

How To: Store and Test Tags



Tags must have sufficient battery power by the MyChild system.

How

To prevent premature battery drain, M (FZB-050 available on our online store [needed-for-tag-storage-sku-fzb-050/](#)). frequency) “chatter” between themselves

Because the tags are always active, to a individual foil pouches should be kept in an RF-shielding solution.

Additionally, the tags should be stored can drain tags’ batteries. Some examples noise are:

- 1) Computers: desktops, laptops, and tablets
- 2) Flat screen and tube-type monitors and

Tag Battery Testing

Must be done for all MyChild

NOTE: Before using the pocket tag reader, the pocket tag reader/tester is powered with the battery compartment. The battery compartment is on the back of the device.

To turn the pocket tester on, press and hold the power button for 4 to 5 seconds.

Once the device is turned on, follow the prompts on the screen.

Step 1

When you see *User Menu* – press the *U* button.

Step 2

The prompt will read *Tag Test Mode*.

Step 3

Press the *ON/ACCEPT* button.

Cut-Band Tag Tamper Testing

Required for Cut-Band Tags Only

Cut-band tags require a tamper test that uses a test plate (CBP-100 available on our online store: [plate-for-cut-band-tag-sku-cbp-100/](https://www.mychildcare.com/plate-for-cut-band-tag-sku-cbp-100/)). If you fail the test, you will time out and have to start again.

Step 1

Place tag on cut-band test plate. Do not touch the tag pins.



After testing a tag, you can assign the tag to a patient at your workstation.

TLM Testing

The pocket tag tester can also be used to test and change the TLM mode.

To turn the pocket tester on, press and hold ON/ACCEPT button **and** the UP button simultaneously for 4 to 5 seconds. Once the device is turned on, follow these steps:

Step 1

You will see Tech Menu. Press the *DOWN* arrow one time.

**If you see the User Menu. Turn the device off and try again, making sure you are pressing both the up button and the ON button.*

Step 2

The prompt will read *TLM Mode*.

Step 3

Press the *ON/ACCEPT* button. This will bring you to the TLM Menu.

Step 4

Press the UP or DOWN button until TLM ENABLE is displayed.

Step 5

Press the Accept button once and RATE MENU will be displayed.

Step 6

Press the Down button until TLM RATE 16 SEC is displayed.

Step 7

Press Accept button.

Step 8

The tag tester displays *Querying TAG*.

Step 9

Hold the tag next to the tag tester to change the TLM mode to be “enabled.”

If successful, the HH will display the Tag’s ID and TLM Enabled

If multiple tags need to be turned on, pressing the Accept button will scan the next tag held close to it.

Make sure there is only one tag at a time near the pocket tag tester while turning on TLM.

Contact Us

For Technical support

TechSupport@McRobertsTech.com

800-776-7328 Option 2

[Download How to Store and Test Tags](#)

For customer service

CustomerService@McRobertsTech.com

800-776-7328 Option 8

Unable to Assign a Tag to a Patient - Troubleshooting



Troub

Topic: Unable to A

Note: Tags auto-enroll when taken out of their RF tags. You must select an available tag as part of the a tag being assigned. If you do not see the tag you

Step 1: After clicking “Enroll Patient” under the drop-down list.

- The tags are always active. Tags should be in pouches when not in use to block their signal. Once removed from a pouch, door controllers will receive the tag’s signal. It’s important to give the system time to recognize the tag once it’s outside the pouch. This can take up to a minute.
- Once seen by the system, that tag should

Available: The status "Available" is assigned to all tags when they are first seen by the system. It indicates that the tag was already enrolled. The status "Available" remains for a time period that is configurable.



What's Next?

Check the tag's last-seen date and time. That information is included in the Tag Inventory.

If the time is not within the last minute, then the system is not seeing the tag.

Proceed to Steps 4 and 5 to determine the reason that the system is not seeing the tag.

Step 4: Test Tag Battery with Pocket Tag Tester

To turn the pocket tester on, press and hold ON/AC these steps:

Step 5: Check to see if Tag Locating Message

The tag must have TLM enabled to be seen by the tester on, press and hold ON/ACCEPT button for 4

Step 1

Once on, you will see Tech Menu displayed on the

Note: If you see the User Menu, turn the device off and press the ON button at the same time.

Step 2

The prompt will read TLM Mode.

Step 3

Press the ON/ACCEPT button. This will bring you to

Select TLM Query to see if TLM is enabled.

If TLM is not enabled



[Click here to download.pdf](#)

How to: Store Radio Frequency Tags



McRoberts
TECHNOLOGIES

How to Store F MyChild, StaffWa

Storing tags

For maximum battery life, RF tags should be stored in foil pouches. Foil pouches experience RF (radio frequency) shielding, which prevents the tags from being active.

To order foil pouches, email [Fulfillment](#)

The tags are always active. To avoid the tags from being active when you do not want them to be, the tags should be stored in a foil pouch. In a non-foil environment, the pouches should be stored in a dry, cool environment.

To order a tag storage box, email [Fulfill](#)

[Download the PDF here.](#)

Tag Management

Tag Warranty/Buy-Back



and

1. Tag Ownership

Tags become the property of Customer w
shipment.

2. Tag Warranty

Tags are guaranteed to function for one y
from date of shipment, the Customer wil
*assumption that customers are storing ta
pouches provided by McRoberts and that
McRoberts.*

3. Tag Buy-Back Program

[Download PDF here!](#)

General Functions

This chapter should cover general user functions

How To: Change your Password

To access your user profile, click on your email address or
On this page you can change your password.

**MyChild**

Home

Patients

Tags

Floorplans ▾

Reports ▾

Admin ▾

Help ▾

Personal Details

First Name:

MyChild6

Last Name:

Administrator

Auto-Logout Delay:

0 seconds ▴ ▾

Permission Level:

SuperAdmin

Role:

Nurse - Admin

Pin:

4578

Save Details

Alarm Notifications

☐ Subscribe to System Alarms

System Alarm Selection

Select System Alarms

☐ Subscribe to Patient Alarms

Patient Alarm Selection

Select Patient Alarms

☐ Subscribe to Device Alarms

Device Alarm Selection

Select Device Alarms

Save Alarm Notification Settings

© 2023 - McRoberts Technologies | b



If you do not know your password, click the reset password button on the login page.

Login

Work Email:

Password:

Log in

[Reset password](#)

How To: Enable 2FA

How To: Update SMS settings

Alarms

Alarm Overview

There are three categories of alarms; Patient, Device, and System. The alarm level indicates the severity of the alarm.

Level 1 alarms notify you about events that could lead to a security breach. Level 1 alarms are indicated by pop-up notifications that disappear after 30 seconds unless you click “View Alarm.” The alarm is displayed in the floorplan view. The floorplan view has priority over all other pages. Make sure you complete whatever you are working on before clicking on “View Alarm” because your work will not be saved.

Level 2 alarms also notify you about events that could lead to a security breach but are more serious than the events that trigger a level 1 alarm. It is recommended that security protocols be followed. There are no pop-up notifications. The floorplan view with the alarm automatically takes over your screen; therefore, whatever you are working on will not be saved. Audio is included.

Level 3 alarms notify you about events that could be security breaches. It is recommended that security protocols be followed. There are no pop-up notifications. The floorplan view with the alarm automatically takes over your screen; therefore, whatever you are working on will not be saved. Audio is included.


[Click here to download the MyChild6 Alarm Chart](#)

Patient Alarms

Patient Alarms are alarms that are triggered by assigned tags under certain conditions. Patient alarms show on your floorplan and point to the device that triggered that alarm or saw the tag last. The alarm includes the patient's name, tag number, and last known location.

Level One Patient Alarms:

Level 1 alarms notify you about events that could lead to a security breach. Level 1 alarms are indicated by pop-up notifications that disappear after 30 seconds unless you click "View Alarm." The alarm is displayed in the floorplan view. The floorplan view has priority over all other pages. Make sure you complete whatever you are working on before clicking on "View Alarm" because your work will not be saved.



TAG NEAR CLOSED DOOR

TAG CHECK-IN MISSED

Level Two Patient Alarms

Level 2 alarms also notify you about events that could lead to a security breach but are more serious than the events that trigger a level 1 alarm. It is recommended that security protocols be followed. There are no pop-up notifications. The floorplan view with the alarm automatically takes over your screen; therefore, whatever you are working on will not be saved. Audio is included.

TAG LOITERING NEAR DOOR: This alarm occurs when an assigned tag that is NOT in Bypass Mode is in a door field for longer than one minute and the door is closed. It replaces the TAG NEAR CLOSED DOOR alarm.

LOW BATTERY TAG ON PATIENT: This alarm occurs when the system receives a low battery message from a tag that is assigned to a patient.

MyChild

Home Patients Tags Floorplans ▼ Reports ▼ Admin ▼ Help ▼ Logout

LABOR AND DELIVERY

McRoberts HOSPITAL

LOW BATTERY TAG ON PATIENT

Change Tag

Device Name: East Hall Nurses Stati

Patient Name: Olt, Alex

Assigned Room:

Tag Number: F1E753

Time of Alarm: 03/28/2022 14:03:57

ACCEPT

tag changed

First, click "ACCEPT" button.

1

Labor and Delivery


© 2022 - McRoberts Technologies | build version: 8122

PATIENT UNAUTHORIZED BY MYCHILD

Level Three Patient Alarms:

TAG TAMPER

TAG NEAR OPEN DOOR



TAG NEAR ELEVATORS DOORS CLOSED

TAG NEAR ELEVATORS DOOR OPEN

TAG IN ELEVATOR

TAG OUT OF PROTECTED AREA

Device Alarms

Device Alarms are alarms that are triggered by door controllers, receivers, elevator controllers, or digital Input/Output devices (I/Os) under certain conditions. Device alarms show on your floorplan and point to the device that triggered the alarm.

System Alarms

System Alarms appear to those with the relevant credentials in their profiles, such as IT users. Other users can gain access through profile customization. System alarms also generate SMART messages to McRoberts technical support so that McRoberts technicians can initiate troubleshooting immediately.

Alarm Chart



Alarms in MyChild6

Patient Alarms are alarms that are triggered by assigned tags under certain conditions. These alarms will pop-up on your floorplan and point to the device that triggered that alarm or saw the patient's name, tag number, and last known location.

Level 1 alarms notify you about events that could lead to a security breach. These alarms will pop-up notifications that disappear after 30 seconds unless you click "View Alarm" in the floorplan view. The floorplan view has priority over all other pages in the system; therefore, whatever you are working on before clicking on "View Alarm" because you will be taken to the alarm details page.

Level 2 alarms also notify you about events that could lead to a security breach. These alarms will pop-up notifications that disappear after 30 seconds unless you click "View Alarm" in the floorplan view. The floorplan view has priority over all other pages in the system; therefore, whatever you are working on will not be saved. Audio is included for these alarms.

Level 3 alarms notify you about events that could be security breaches. These alarms will pop-up notifications that disappear after 30 seconds unless you click "View Alarm" in the floorplan view. The floorplan view has priority over all other pages in the system; therefore, whatever you are working on will not be saved. Audio is included for these alarms.

Alarm Name	Level	Alarm Description	Setting Configuration
TAG NEAR CLOSED DOOR	1	This alarm occurs when an assigned tag that is NOT in Bypass Mode is in a protected door field and the door is closed. This alarm self-accepts when the tag is no longer in the door field. If the door opens while the tag is in the field, this alarm gets replaced by a TAG NEAR OPEN DOOR alarm. If the door remains closed but the tag is still in the field for more than one minute, this alarm is replaced by a TAG LOITERING NEAR DOOR alarm.	Yes. You can enable or disable this alarm at individual doors.



Alarm Name	Level	Alarm Description	Setting Config
TAG LOITERING NEAR DOOR	2	This alarm occurs when an assigned tag that is NOT in Bypass Mode is in a door field for longer than one minute and the door is closed. It replaces the TAG NEAR CLOSED DOOR alarm.	Yes. Y disabl indivi doors
TAG NEAR OPEN DOOR	3	This alarm occurs when an assigned tag that is NOT in Bypass Mode is in the field of an open door.	No
TAG NEAR ELEVATORS - DOORS CLOSED	3	This alarm occurs when an assigned tag NOT in Bypass Mode is in a bank style elevator field and the elevator doors are all closed. In a standard design, the elevators do not stop on the floor where the tag is when this alarm is active.	No
TAG NEAR ELEVATORS - DOOR OPEN	3	This alarm replaces "TAG NEAR ELEVATORS - DOORS CLOSED." This alarm occurs when an assigned tag NOT in Bypass Mode is near bank style elevators, and a motion detector at an elevator door detects motion, thereby indicating the elevator door is open. In a standard design, the elevators are disabled on the floor where the tag is when this alarm is active.	No
TAG IN ELEVATOR	3	This alarm occurs when an assigned tag NOT in Bypass Mode is in the field of an individual elevator controller. In a standard design, the elevator door will stay open, and the elevator will not move when this alarm is active.	No



Alarm Name	Level	Alarm Description	Setting Configurable
PATIENT EXIT UNAUTHORIZED BY MYCHILD	2	This alarm occurs when someone enters a correct PIN code at a door and there is an assigned tag in the door field and that tag is NOT in Bypass Mode. The egress is authorized by the hospital because the correct PIN was entered; however, the egress is unauthorized by MyChild.	This alarm can be turned on or off for individual doors.
TAG OUTSIDE OF PROTECTED AREA	3	This alarm occurs when a tag is picked up by a receiver located outside of a protected unit/zone.	No
TAG TAMPER	3	This alarm occurs when the circuit formed by the connection of a cut-band and a cut-band tag is interrupted OR the sensing mechanism on the tag no longer senses that it is on skin.	No
TAG SIGNAL MISSED	1	This alarm is triggered when the system does not receive any TLM signals from an assigned tag for a configurable time period.	Yes. You can choose how many seconds should elapse without hearing from an assigned tag to trigger this alarm. The default is 30 seconds.
TAG SIGNAL MISSED 2	2	This alarm is raised after a configurable period of time after a Tag Signal Missed.	Yes. You can choose how many seconds should elapse after a Tag Signal Missed.



Bypass Alarms: MyChild provides a ByPass Mode that allows you to take a patient without setting off alarms. You must select a departing unit and a destination. When a patient status changes from “Active” to “Bypass” when a patient is in ByPass Mode.

When the time allotted for the bypass is almost expired, a “BYPASS EXPIRING” alarm is triggered. When the timer reaches zero, a “BYPASS EXPIRED” alarm is triggered. All bypass alarms can be edited by you to edit the bypass. Changes can also be made on the patient card. When the alarm is “Complete,” you are acknowledging that the patient has arrived at the destination.

Alarm Name	Level	Alarm Description	Setting Configurable
BYPASS STARTED	1	The first bypass-designated door passed by a patient in bypass triggers this alarm and starts the bypass timer.	No
BYPASS EXPIRING	1	This alarm occurs five minutes before a patient in bypass is due to pass the second bypass door (destination). If the patient passes the second bypass door during the five-minute period, this alarm is replaced by a “BYPASS COMPLETE” alarm. If five minutes pass and the patient has not passed the second door, this alarm is replaced by a “BYPASS EXPIRED” alarm.	You may select how many minutes before the bypass timer expires to raise this alarm. The default is five minutes.
BYPASS COMPLETE	1	This alarm occurs when the patient in bypass passes the second bypass door.	No



Suspend Tamper Alarms: MyChild provides a suspend tamper function v or sensing tag temporarily without setting off alarms. When you put a pa patient status changes from “Active” to “Tamper Suspended.”

When the time allotted for the suspend tamper is almost expired, a “SUS triggered. When the timer reaches zero, a “SUSPEND TAMPER EXPIRED” alarms allow you to edit the suspend tamper. Changes can also be made

When you click “Complete” to complete a Suspend Tamper, you are ackn patient.

Alarm Name	Level	Alarm Description	Setting Configurab
SUSPEND TAMPER STARTED	1	This alarm notifies you that a patient has been put into Suspend Tamper Mode and that the Suspend Tamper timer has started. The patient status changes from “Active” to “Tamper Suspended.”	No
SUSPEND TAMPER EXPIRING	1	This alarm appears five minutes before the suspend tamper allotted time is due to expire. At this point you can edit the tamper suspend to add time or complete the tamper suspend.	You may select how many minutes before the timer reaches ; to raise this alarm.



Device Alarms are alarms that are triggered by door controllers, receiver Input/Output devices (I/Os) under certain conditions. Device alarms show the device that triggered the alarm.

Level 1 alarms notify you about events that could lead to a security breach. Level 1 alarms show pop-up notifications that disappear after 30 seconds unless you click on them. Level 1 alarms are displayed in the floorplan view. Floorplan view has priority over all other views. Therefore, whatever you are working on before clicking on a level 1 alarm because the floorplan view takes over your screen.

Level 2 alarms also notify you about events that could lead to a security breach. Level 2 alarms are the events that trigger a level 1 alarm. There are no pop-up notifications for level 2 alarms. A level 2 alarm automatically takes over your screen; therefore, whatever you are working on before clicking on a level 2 alarm because the floorplan view takes over your screen. Audio is included.

Level 3 alarms notify you about events that could be security breaches. Level 3 alarms are the events that trigger a level 2 alarm. There are no pop-up notifications for level 3 alarms. A level 3 alarm automatically takes over your screen; therefore, whatever you are working on will not be visible. Audio is included.

Alarm Name	Level	Alarm Description	Setting Configurable
LOW BATTERY TAG IN INVENTORY	1	This alarm is triggered when a tag that is not assigned to a patient sends a Tag Initiated Communication (TIC) message of low battery. The closest device to this tag reports this.	You can change the number of TIC messages received before triggering this alarm.
GLOBAL LOCKDOWN OFFLINE	2	This is for I/O devices.	The number of signals missed within a certain time period can be configured.



Alarm Name	Level	Alarm Description	Setting C
DOOR CONTROLLER OFFLINE	2	MyChild continuously requests status from each door controller. No response from a door controller triggers this alarm. When a door controller is offline, it still locks when a tag is in the field; however, users will not see door-related alarms on their workstations. Of course, if a door is offline due to a loss of power it will not lock when a tag is in the field.	You can c number c before tri alarm. Th signals (a seconds).
BYPASS DEVICE OFFLINE	2	This is for doors that are designated Bypass Doors and tag detectors utilized for Bypass.	You can c number c before tri alarm. Th signals (a seconds).
INVALID PIN ENTERED	2	This alarm occurs when too many incorrect bypass PIN codes are entered on a MyChild keypad.	Yes. The attempts which the can be m configure 3 tries in
ANTENNA DISCONNECTED	3	This alarm means that doors will not lock! MyChild software continuously requests status from each door controller. If an antenna is disconnected, this message is included with the status, triggering	You can s of exciter disconne within a c period to alarm. Th



Alarm Name	Level	Alarm Description	Setting Configuration
WIRELESS RECEIVER OFFLINE	1	MyChild software continuously requests status from each device. No response from a wireless receiver triggers this alarm. This is usually a temporary network issue.	Yes. You can select how many signals are missed in a configurable time period. The default is 30 signals in 30 seconds.
RECEIVER OFFLINE	1	MyChild software continuously requests status from each device. No response from a wired receiver triggers this alarm. This is usually a temporary network issue.	Yes. You can select how many signals are missed in a configurable time period. The default is 30 signals in 30 seconds.



System Alarms appear to those with the relevant credentials in their p can gain access through profile customization. System alarms also gen technical support so that McRoberts technicians can initiate troublesho

Alarm Name	Level	Alarm Description	Setting Configurable?
NETWORK IN DEGRADED STATE	1	This alarm occurs when, during a five-minute period, there are three or more ping failures OR average ping latency is 0.5 milliseconds or greater. A ping is an echo request sent from the server to the INM. Latency is the amount of time it takes for the request to occur measured in milliseconds.	No
NETWORK OFFLINE	2	This alarm occurs when there is an interruption in power or communications in the MyChild network. This is a security risk because without connection to the server, MyChild loses its functionality, so users don't see important alarms. This condition does not impact physical security. The doors operate independently and still lock and alarm when a tag is in proximity.	No



Level 0 alarms are notifications that do not indicate a security risk and

Alarm Name	Level	Alarm Description
TAG AUTO-ENROLLED IN SYSTEM	0	This notification is logged when a tag auto-en
PATIENT ENROLLED	0	This notification is logged when a patient enr
PATIENT UPDATED	0	This notification is logged when patient infor
PATIENT UNENROLLED	0	This notification is logged when a patient is u
PATIENT LOCATION DISPLAYED	0	This notification is logged when the LOCATE f
BYPASS CREATED	0	This notification is logged when a bypass is cr start the bypass timer; however, the patient :
DOOR BYPASSED	0	This message occurs as part of the normal bypa bypass door. This initiates the bypass timer.
PATIENT BYPASS ENDED	0	This notification is logged when a patient pas completes the bypass.
SUSPEND TAMPER COMPLETED	0	This notification is logged when the suspend status goes back to "Active."
DEVICE BUFFER OVERFLOW	0	This notification is triggered by the devices whe This is self-healing and no user action is require
PIN ENTRY ACCEPTED	0	This notification is triggered anytime a valid PIN this is logged as an informational message.
DEVICE PIN CHECKSUM FAILURE	0	This notification is triggered when the door con internal PIN database. This is normally a self-he where power cycling the device and/or restarti
CONTROLLER /	0	These notifications are normally triggered imm

[Click here to download the Alarms Chart](#)

Spec Sheets

Cut-Band Spec Sheet



McRoberts
TECHNOLOGIES

SPE
PRO
SKU

APPLICATIONS USED IN:

- MyChild
- WanderWatch

OVERVIEW:

The Cut-Band is a security and monitoring device consisting of a thin, soft and flexible band of conductive material. The tag has metal pins on either side that create a loop that can be pulled through. The bands are one-way adjustable and can be removed by being cut off.

FEATURES:

- Tag Locating Message (TLM) provides tag status
- If TLM signals are missed, an alarm is generated
- Automatically engages doors and elevators
- Cutting the band to remove a tag generates an alarm
- Option for sensing capability is available.

[Download the Cut-Band SpecSheet](#)

Tummy Tag Spec Sheet



McRoberts

TECHNOLOGIES

APPLICATIONS USED IN:

- MyChild

FEATURES:

- Tag Locating Message (TLM) provide
- If TLM signals are missed, an alarm i
- Tamper deterrent umbilical attachm
- Lightweight and compact design.

[Download the Tummy Tag SpecSheet](#)

Bracelet Tag Spec Sheet



McRoberts
TECHNOLOGIES

SPEC S

SKU: B

APPLICATIONS USED IN:

- MyChild
- WanderWatch

OVERVIEW:

The Bracelet Tag is a security and monitoring device is applied to patients using a bracelet. The small, light is applied to the patient using a variety of commonly a

FEATURES:

- Tag Locating Message (TLM) provides tag su
- If TLM signals are missed, an alarm is genera
- Automatically engages doors and elevators f
- Lightweight and compact design.
- Unique electronic serial number identifies ea
- Tags auto-enroll in the system.
- High-capacity 3V lithium ion coin cell batter
- Eastman Tritan TX1001 amorphous copolyes
- Smooth surface for easy cleaning

[Download the Bracelet Tag SpecSheet](#)

Wireless Receiver Spec Sheet



McRoberts

TECHNOLOGIES

APPLICATIONS USED IN:

- MyChild
- AssetWatch
- StaffWatch
- WanderWatch

FEATURES:

- Stand-alone devices receive tag signal through standard Wi-Fi infrastructure.
- Powered by three-pronged wall socket.
- MAC address enables customer to identify devices.
- Receiver coverage customizable to customer's needs if desired.

[Download the Wireless Receiver SpecSheet](#)

Bypass Keypad Spec Sheet



McRoberts
TECHNOLOGIES

SPEC
SKU: ,

APPLICATIONS USED IN:

- MyChild
- StaffWatch
- WanderWatch
- UrgenTag

OVERVIEW:

The Bypass Keypad is installed close to an exit equipment controller. A staff member can, by entering his/her ID into the door controller field with a tag without generating an alarm event and the identity of the staff member is recorded.

FEATURES:

- Accommodates up to 1000 individual user Profiles
- 3 LED indicators: Standby, Alarm, Bypass.
- Audible alarm (adjustable volume).
- Operates even if communication is lost with the controller.

[Download the ByPass KeyPad SpecSheet](#)

STEP Customers - Instructions



Simple Tag Exchange

At the STEP program initiation, McRoberts contract. Included in the package will be box) and a return shipping label. You will McRoberts.

Each tag will arrive in its own foil pouch. foil pouches in a metal tag storage box so McRoberts Storage Box can accommodate

Each time a tag is removed from a patient weekly basis, send the used-tag box to Mc tag box inside a shipping box. Put the ret

Ship the box to:

McRoberts Technologies
STEP Fulfillment
Suite 111
409 Joyce Kilmer Avenue

[Click Here to Download the PDF](#)

Changelogs

This chapter holds changelogs for each newly released version.

v6.0.7-SP2

Release notes - MyChild6 - v6.0.7-SP2

Task

MST6-497 Login and Authentication - Group Permissions

Description

The Problem: There is currently no way for end-users or system administrators to easily update user roles and permissions.

Recommended Resolution: Create a new form/set of forms that allow the systems administrator to configure the system permissions.

MST6-801 Enable Windows Auth for SQL Server

Description

Problem: While deploying A customer we discovered that they only have Windows Auth enabled for their central SQL Server. We have requested to have SQL Server Auth enabled, but after doing some research we found that enabling Windows auth really should not be difficult. Below is documentation of what changes need to be made. We want to enable this so that in the future we can use Windows auth / LDAP more easily.

Documentation:

- [Configure Windows Authentication in ASP.NET Core](#)

MST6-809 Patient functions Reorder

Description

Our patient card has 4 buttons on the back, also known as functions. We need to reorganize them so they have a hierarchy.

The order of the buttons will display in this order:

1. BYPASS
2. SUSPEND TAMPER

- 3. LOCATE
- 4. UNENROLL

EDIT: additional corrections

- Live map broken overflows have been removed.
- Filter on directory view corrected
- width of device config page corrected.
- tabstrip-1/5 adjusted with selector call []

MST6-843 Database needs more indexing.

Description

AS A: Developer

I NEED: to optimize the backend database queries and indexes

SO THAT: data can be retrieved from the database faster and more efficiently

MST6-847 Create a new setting for Tag Signal Missed Delay Period

Description

AS A: developer OR system admin

I WOULD LIKE: to add a new setting to the FacilitySettings that will allow the facility to control the amount of time between each Tag Signal Missed Alarm

SO THAT: the facility can lessen the number of TSM alarms raised, assuming they have lots of tags and patients

Further Discussion:

We have found that at larger facilities, the number of TSM alarms raised per hour can get to be high, borderline uncontrollable, this is mostly because we repeat the alarm every 30 minutes until we start hearing from the tag again.... By adding a new setting, we could give the facility the ability to delay the repeated alarms for XX number of minutes (maximum of 9999)

MST6-859 Improve Device & Tag Caching for Wireless Plug Messages

Description

Problem: We are still not processing events as fast/efficiently as we need to be.

Story

MST6-576 Account - system lock after several failed login attempts

Description

AS A user

I NEED the system to lock after several failed attempts

SO THAT I can reset and account. **OR** prevent unauthorized people from using the software.

The question we ran into was this:

"Does the system lock-out an ID after several failed attempts at login? If yes, how is the ID unlocked?"

*More information here ****

MST6-613 Advanced Filter/Sorting Styling

Description

Objective: Now that we have our filter and sorting working as we would like, it's time to style the filter.

Currently, the filter rests in a navbar above the card view on the far-left side. As we show progress on this page, we may run into some adjustments from the team.

DEV EDIT:

removed the ridge style a while ago which made this look better on all solutions. No need to move the position from the top left unless requested by customers.

MST6-779 Dashboard | Unaccepted Alarms | Scrollbar and Polish

Description

We feel that the unaccepted alarms tile on our dashboard can be polished with styling to make it look and feel better. When resizing the window to be small or on mobile applications, text compresses and does not look too good, specifically in the tabs. We also feel the scrollbars on this tile are clunky and unattractive.

OBJECTIVE: Adjust the scroll bar to be more visually appealing. This can be done with color adjustments as well as size adjustments. We must verify that adjustments to any scroll bar work ACROSS ALL BROWERS INCLUDING IOS AND ANDROID systems. Bootstrap also has scroll bar adjustments we can use.

We also want to prevent tabs and text from wrapping to new lines when windows get resized.

MST6-781 User Management | User Control | input field & Active button

Description

12/28/23 During our mychild meeting we found some things we would like to change on our User Management page.

First, with our User Control popup where we make an account, we need to correct the PIN field to

match all the others. It needs to be taller and grey.

Second, the active checkbox must be set to active initially. It is a safe assumption that they will be used if a user makes accounts.

On the main view, (Grid) we have commands or buttons per account. One of the button's texts needs to be adjusted. Our disabled account should act as a toggle, with the text of "activate" and "Deactivate". this keeps actions and text consistent.

MST6-782 EnrollPatient | CreateBypass | Field Reorder & Refactor

Description

It was requested that reorganize our enroll patient popup. We started a discussion to reorder some things to make more sense, an example was to move the sex field up by the patient details. This discussion turned into a pool of ideas, but it was decided that we would like two sections to this popup. One for location details such as tag and device, and one for patient information.

We also would like to see a better / bigger enroll button on this popup. The team would like to see some shaking animations or hover effects on the button.

I need to create mocks to provide visuals to continue this discussion. It was requested to make this popup more appealing overall. When mocks are made and approved, we can start development.

Similar story item:

"Our Patient Bypass popup currently works but the team feels it's slightly disorganized. We have patient info input fields at the top and the bottom, it just felt messy.

I brought an idea to rebuild the Enroll popup, however the team felt it was overkill as it really does not need to be rebuilt.

WE MUST re-organize the fields on the popup. All patient information will be moved to the top (including the sex checkbox). Tag information and location will be at the bottom of the popup.

The string of Department must be changed to Unit.

The create bypass button will need to change to "Put Patient in Bypass". "

MST6-806 Correct Device Config page for ALL solutions

Description

In the merge project, each solution has a broken device config page when loading in a map/floorplan for the first time. this should be corrected using a template.

This is a good starting point for templates for each solution using Brad's environment variable.

MST6-818 Simulator UI & Functionality broken.

Description

simulator.dev.(URL).com is showing some crazy UI. Large carrots and broken cards.

MST6-835 Viewing and Exporting reports is slow if not impossible.

Description

The Problem: When attempting to view or export reports, almost all of them will either fail to load, or fail to export, and doing either takes quite a long time, and the longer the system runs the longer these reports take to run.

Recommended Resolution: review the report implementation and refactor the backend queries to retrieve a minimal amount of data needed for the reports.

MST6-840 Map views and corrections

Description

Maps and the related pages have some style errors that need to be resolved. Now that the projects are merged, we need to ensure that the UI projects display the map elements correctly.

MST6-852 Trade Show | Demo Bugs | Inconsistencies

Description

During and after the trade show in Orlando, we found several small adjustments we would like to make to our StaffWatch solution, along with a handful of minor corrections.

- Map upload and device configuration have a handful of minor issues such as padding, color, spacing, overflow, and text that need to be corrected
- The password reset page needs to be rebuilt
- The logout page should say the solution name, not only "Mychild"
- Headers are not uniform on the directory card view for all solutions. Mychild seems to be perfect, the others do not seem right and look off.
- The close button on most popups not showing correctly. Telerik error.
- Assetwatch assign button is not fully clickable due to the logo container not hiding when we flip the card.
- StaffWatch has "patient" alarm options on the MyProfile page.
- Logos not working on MyProfile page.

MST6-858 No Ability to edit/delete User Groups

Description

AS A: user OR system admin OR user manager

I NEED: the ability to add/edit/delete user groups, we currently can ADD new user groups, but we cannot edit or delete the groups once they have been created

SO THAT: we can further customize the software to meet our individual needs

Bug

MST6-820 Alarm Times | Wrong day displaying.

Description

It was found that the alarms were displayed at the correct time but the incorrect day.

MST6-821 Unassigned Tag | TIC alarms

Description

Problem: When a cut-band tag is unassigned, Tamper alarms are being raised

Recommended Resolution: When a Cut-Band or Tamper tag is *NOT* enrolled, then TIC alarms should be ignored as the tag is more than likely in storage.

MST6-826 Level 3 alarm audio being intermittent.

Description

It was found the alarm sound for level 3 does not always go off. We need to find a solution for this bug.

-Correction for this issue will be attached to MST6-820

MST6-827 Duplicate Serial Numbers in Database

Description

Duplicate serial numbers should **NOT** be allowed in the database.

MST6-828 Deleting Devices

Description

Deleting a device fails if the device icon has not been removed from the map first.

MST6-830 Shared MenuLayout throwing InnerHTML error.

Description

Problem: when viewing the main dashboard, a JavaScript error is being logged that is saying it cannot set “inner HTML” in function “poll”

Recommended Resolution: review the poll function and verify that the function is no longer used and remove it.

MST6-832 System struggles when system has more than 100 devices & tags.

Description

The Problem: We have found that when a system has > 100 devices and/or >100 tags, the system seems to become unresponsive, especially if the database server is shared between all 3 applications, if one of the 3 applications has a high number of devices, then all applications will suffer degraded performance.

Resolution: Review the “TagEventReceived” code and optimize the code to use caching & smaller, quicker SQL queries to attempt to retrieve data more efficiently and to help improve overall system performance

MST6-833 Small UI Corrections | Hutch

Description

We have a handful of corrections across the watch solutions.

ASSETWATCH:

- Edit asset now has an 'X' to close the popup instead of the broken '?'
- Broken border 'Ridge' style on tag page removed.

ALL WATCH SOLUTIONS:

- Alarm popup no longer has an 'x' floating on the righthand side.
- Tag Locate All page is now dynamic to the solution being used. That included text correction in the tab view of the grid.
- Mychild will display 'Patients' and headers will be turquoise
- Assetwatch will display 'Assets' and headers will be dark blue
- Staffwatch will display “Staff” and headers will be pink
- Wanderwatch will display “Resident” and headers will be blue

plus a handful of minor corrections.

MST6-834 Data Queries are pulling far too much data.

Description

The Problem: While investigating performance issues we noticed that some of our database queries are joining up to 20 tables and selecting up to 100 columns of data at once, additionally

some of the queries are requesting the entire table worth of data, which for some tables can be upwards of a million rows, or more.

Recommended Resolution: review the database queries and refine them so that only the required data is pulled, attempt to minimize the number of columns and rows queried for all queries, there is almost never a reason to pull the entire table worth of data, especially when you are talking about a log table that could be millions of rows in size.

MST6-836 Tags Not being removed from the system.

Description

The Problem: For some reason, we are not only seeing quite a few non-McTech tags being added into the system, but we are also having an issue where tags can be removed from the physical area, but they are never removed from the McTech system, even after being “missing” for months, they still appear as active in the system

Recommended Resolution: Review the AutoEnroll and AutoUnEnroll code to ensure that tags are disabled and removed on the correct schedule.

MST6-837 AutoUnEnrollService is still ignoring AutoEnrollPending and Disabled tags.

Description

The Problem: The AutoUnenroll service is still ignoring the primary group of tags that should be removed.

Resolution: the base SQL query for the unenroll service was inadvertently filtering out the AutoEnrollPending and Disabled tags, so the service was never seeing any tags for removal

MST6-838 Clicking alarms on missed alarms causes a JavaScript error.

Description

The Problem: When clicking on an unaccepted alarm in the unaccepted alarms list, a JavaScript error is thrown instead of the alarm being displayed.

Resolution: the underlying dataset query is not returning the correct data to display the alarm, which is valuing the underlying JavaScript to be invalid

MST6-841 Tag List does not allow sorting by LastSeenTime

Description

Problem: The Tag list does not allow sorting by LastSeenDateTime

Resolution: include the LastSeenDateTime in the sort descriptors

MST6-842 McPlugTag Events are all being ignored.

Description

The Problem: MC6/SW6/WW6 are all blocking and discarding almost all McPlug Tag Messages, without processing.

Resolution: After review, we found that the RSSI Filter that was implemented for AssetWatch, was incorrectly implemented and was causing the filter to skip all events unless you were running AssetWatch AND the RSSI value was greater than -80

MST6-844 Tag Signal Missed Alarms are out of control.

Description

The Problem: for some reason, when a patient/asset/staff/resident is enrolled to a tag, the tag will go into nearly non-stop tsm mode raising tons of alarms for what appears to be no reason.

Recommended Resolution: review the Tag Signal Missed generation code to determine why so many alarms are being raised.

MST6-848 Tag Unenroll is still not functioning correctly.

Description

The Problem: Tags that have not ever been seen (IE: LastSeenDateTime is all zeros) are never being removed from the system, which is causing confusion about which tags are usable and which are not.

Recommended Resolution: Review the Unenrollment code and determine why tags are not being unenrolled from the system.

MST6-849 Tag Signal Missed Alarms generated when they should not be.

Description

Problem: We are seeing quite a few Tag Signal Missed alarms being raised for tags that *ARE* checking in correctly and promptly. It seems like only the tags that are correctly checking in are having TSM alarms raised.

Recommended Resolution: Review the TSM Alarm Generation code to determine why TSM alarms are being generated when they should not be.

MST6-850 Tag Caching is only partially working.

Description

The Problem: The new Tag caching algorithms are working well, however, we discovered that there is at least one place where the logic could cause heavy queries against the dataset.

Resolution: Refactor the null logic to create an in-memory cache entry with a status of “Locked”, this should be cleared every minute OR by the auto-enrollment process upon successful autoenrollment of the tag

MST6-851 Reports are showing Numbers or Blanks instead of Patient And/or Tag Names/Numbers

Description

Problem: When viewing the Unfiltered Log Report there are several instances where the Tag Number and Patient Name are blank

Resolution: Review the Reporting code and determine why patient/staff/asset names and tag numbers are not being properly reported

MST6-853 Changes in Drop Downs in MyChild List View Reverting on Page Refresh

Description

issue Description: Testing editing patients in the list view of MyChild6 shows that any drop down changed in list view reverts back to original setting after a page refresh. This does not occur in other RTLS 6 applications. Nor does this effect text entries and check boxes.

Replication Instructions: In list view edit a patient. Change a drop-down selection. Click save. Refresh page.

Functionality Suggestion: No setting should ever revert to page refresh.

Browser: Chrome

OS: MacOS

Build: Dev

MST6-854 Incorrect Item Count on Tags Page in MyChild

Description

Issue Description: The number of item counts at the bottom of the Tags page in MyChild is not showing the correct number of tags listed on the page.

Replication Instructions: N/A

Functionality Suggestion: The number of items should match the number of tags on the page.

Browser: Chrome

OS: MacOS

Build: Dev

MST6-855 Database Updates causing delays and waits.

Description

The Problem: With a large system, with more than 200+ devices and tags, we are seeing an issue where updating the LSDT on the Tags and the LH on the Devices is causing very heavy loads on the SQL server, often causing the CPU usage to be upwards of 60%.

Resolution: Review the underlying database management code and determine if any improvements can be made to the update process

MST6-856 No ability to edit/delete User Group Permissions

Description

Problem: if you browse to the sysadmin section, you can go to the Identity section which gives you the ability to add new groups, as well as add permissions to the groups, however, there is currently no way to edit or delete permissions once they have been assigned to the group.

Resolution: Review the Group Permission Manager under sysadmin and determine why edit/delete are not options

MST6-857 Changes to Patient Floor or Department not saved in ListView

Description

issue Description: Testing editing patients in the list view of MyChild6 shows that any drop-down change in list view reverts to the original setting after a page refresh. This does not occur in other RTLS 6 applications. Nor does this affect text entries and checkboxes.

Replication Instructions: In list view edit a patient. Change the drop-down selection. Click save. Refresh page.

Functionality Suggestion: No setting should ever revert to page refresh.

Browser: Chrome

OS: MacOS

Build: Dev

6.0.7-SP1

Release notes - MyChild6 - v6.0.7-SP1

Bug

MST6-810 DataProtection Keys are not loaded properly

Description

The Problem: during the merge project it looks like the DataProtection configuration got changed, when upgrading an existing Windows install, the patient names are all encrypted and you have to add the DataProtection configuration to fix it.

Recommended Resolution: Correct the data protection loaders so that the keys are correctly loaded during startup

MST6-811 LDAP users cannot log in with user-friendly domain

Description

The Problem: When configuring LDAP integration, we have to configure the application to query the AD server directly using the full domain name (IE: [ad01.mydomain.com](#)), however, users will typically have a username like [myuser@mydomain.com](#). This is a problem because the code expects you to log in with the same “domain” that is set in the appsettings.json

Recommended Resolution: update the backend logic to check to see if the main parts of the domain name match and then process the request based on that result (IE: [user@mydomain.com](#) should be the same as [user@ad01.mydomain.com](#))

MST6-812 Sys-Admin | Alarm Acceptor

Description

It was found our mass alarm acceptor in our system admin page is not working as it should.

I have found a 404 error and a reference error on the index page.

MST6-813 LDAP User roles are not updated

Description

The Problem: When a users group (IE: Role) is changed in LDAP, that change is not reflected in the MyChild Software, so for example, if a user is created with a role of “ClinicalStaff” and then at some point down the road they are placed into a different group (IE: ClinicalAdmin) the role is never updated in MyChild

Recommended Resolution: When a user logs into the system, we should verify that thier role matches whatever LDAP group they are in, if they don’t match (IE: LDAP is different from MyChild), then we should updated the MyChild side to match the LDAP side

MST6-814 LDAP Disabled Accounts Able to login

Description

The Problem: When an account has been disabled in the LDAP Server, the account should *NOT* be allowed to login to the software anymore, however for some reason the LDAP active flag seems to be being ignored for some reason.

Recommended Resolution: update the login procedure to verify the account active flag, if the account is not active (IE: disabled) then the login should be denied

MST6-815 MyProfile | Static size

Description

Our MyProfile page will break when a user zooms into the page using the default zoom on the web browser (ctrl + ScrollWheel).

Our navigation menu has an overlay that starts to break and push things around.

All pages attached to the myprofile page need to use percentages for dynamic sizing, this included the text fields as they may overlap with Zoom.