



Reporting and Analytic Capabilities of MyChild6

Reporting and analytics are key features of MyChild6. All data is available in “Reports.” All data can be sorted to your preference and/or exported to Excel for analysis. Reports can be set up to be auto-emailed to you at a frequency of your choosing.

Reports are accessed through the main menu under the “Reports” tab.

Reports include:

Tag Reports

- Inventory: Shows all the tags you currently have in your inventory. For each tag you will see the tag number, type, status, and last seen location with a date and time.
- Activity: Select a specific tag from the drop-down list and see all activity for that tag. This includes when the tag was auto-enrolled in the system.
- Alarms: Select a specific tag from the drop-down list and see all alarms that tag has triggered.

Patient Reports

- Directory: Shows all the patients you currently have enrolled in the system. For each patient you will see the tag number, status, and last seen date and time.
- Activity: Select a Patient from the drop-down list and see all activity for that patient. Patient activity that does not trigger an alarm can be found here.
- Alarms: Select a patient from the drop-down list and see all alarms that have been triggered for that patient. This includes information on who accepted the alarm and the annotation.

Device Reports

- Inventory: Shows all the devices you currently have in the system with their serial numbers, types, names, and selected settings.
- Activity: Select a device from the drop-down list to view all activity for that device.
- Alarms: Select a device from the drop-down list to view all alarms triggered by that device.

User Reports

- Directory: Shows all users who have access to your MyChild system.
- Activity: Select a user from the drop-down list to view all activity for that user.

Log

- Clinical: Shows all activity and alarms for anything having to do with a patient or a tag.
- System: Shows all activity and alarms that have to do with the system itself.
- Unfiltered: Shows everything from every type of report.

MyChild6 also includes a feature called System Monitoring and Remote Technology (SMART). SMART automatically emails McRoberts Technologies when system alarms are triggered. Users can also receive these notifications by selecting this feature in their settings. Each user may select email or SMS as the alarm notification communication preference.