

# How To: Enroll a Patient

Welcome, [admin@mychildcloud.com](#) | 10 | [Show](#)

Powered by McRoberts Technologies [Visit us to](#)

Home Patients **Click Patients Tab** Admin Help Logout

### TAG SUMMARY

| Status     | Count |
|------------|-------|
| InUse      | 6     |
| Available  | 44    |
| LowBattery | 21    |
| Bypass     | 1     |

Inventory Target: N/A Inventory Total: 72


### SUMMARY

| Status | Count |
|--------|-------|
| Active | 6     |
| Bypass | 1     |

### UNACCEPTED ALARMS

|  |         |
|--|---------|
| Tag near open door - 02/11/2022 20:11:32 | Level 3 |
| Tag near open door - 02/11/2022 20:09:32 | Level 2 |
| Tag near open door - 02/11/2022 18:44:03 | Level 1 |
| Tag near open door - 02/09/2022 21:04:40 |         |
| Tag near open door - 02/09/2022 21:04:35 |         |
| Tag near open door - 02/09/2022 21:04:31 |         |
| Tag near open door - 02/09/2022 21:04:26 |         |
| Tag near open door - 02/09/2022 21:04:21 |         |
| Tag near open door - 02/09/2022 21:04:18 |         |

[dev.mychildcloud.com:8080/Patients/Directory/](#)



The enrollment window will pop up. Select the tag you wish to use from the drop-down list.

**Enroll Patient** [X]

First Name

Middle Name

Last Name

Room

Gender

☐ Flight Risk

Unit

Tag

Floor

© 2024 - McRoberts Technologies | build version: 1.0.0

Home Patients Tags Floorplans Reports Admin Help My Profile

**Enter the patient's information. If you leave any required fields blank you will be prompted to fill them out. Only the "Middle Name" field is optional and can be left blank.**

**Enter Patient Information**

**Then Click Enroll**

**Enroll Patient**

First Name

Middle Name

Last Name

Room

Unit

Tag

Floor

Sex ☒ Male ☐ Female

☒ Flight Risk

☒ Mother Required

☐ Court Order

☐ Other

© 2024 - McRoberts Technologies | build version: 1.0.0

Patient will now have a status of "Active" and the total number of patients will be displayed on the dashboard. Other patient status are "Bypass" and "Tamper Suspended."

---

Revision #7

Created 1 February 2022 16:54:05 by Admin

Updated 20 February 2025 15:38:40 by Anna Olt