



APPLICATIONS USED IN:

- MyChild
- WanderWatch

OVERVIEW:

The Bracelet Tag is a security and monitoring device that is applied to patients using a bracelet. The small, lightweight tags can be applied to the patient using a variety of commonly available bands.



FEATURES:

- Tag Locating Message (TLM) provides tag supervision and locating.
- If TLM signals are missed, an alarm is generated.
- Automatically engages doors and elevators for perimeter security.
- Lightweight and compact design.
- Unique electronic serial number identifies each tag.
- Tags auto-enroll in the system.
- High-capacity 3V lithium ion coin cell battery. One-year prorated warranty.
- Eastman Tritan TX1001 amorphous copolyester case provides:
 - Smooth surface for easy cleaning
 - Attractive, high-impact case
 - Hydraulic stability
 - Chemical resistance
 - Heat resistance

SPECIFICATIONS:

Operating Temperature:	14°F (-10°C) to 122°F (50°C)
Humidity:	95% non-condensing
Dimensions:	1.1" x 1.0" x 0.5" (2.79cm x 2.54cm x 1.27cm)
Weight:	0.3 oz (8.5g)
Input Frequency:	307kHz
Output Frequency:	433.92MHz
Transmit Rate:	16 Seconds for Supervision. Immediate alarms.